## Citra-Lock™

## Catheter-lock 20 x 5 ml vials

## Composition:

46,7% Trisodium Citrate in water. The Citra-Lock is delivered sterile and supplied as a clear solution. The pH is adjusted to approximately 6,7. Aseptic technique must be maintained at all times. Each single-use vial contains 5 ml.

## Prescription:

Citra-Lock is indicated as a catheter lock solution to prevent coagulation of the blood and infection in any type of intravenous catheter. To prevent:

- Clotting of blood in catheter
- Bacterial infections in catheter

## And reduction of:

- Biofilm formation in catheter

## Performance:

Trisodium Citrate (TSC) causes anticoagulation by chelation of ionized calcium in the extracorporeal circuit into a soluble complex. Because calcium is an integral ion involved in the clotting cascade, local removal by citrate prevents the activation of clotting cofactors, factor X, and prothrombin, and the ultimate formation of fibrin.

Antimicrobial effect and a reduction of biofilm formation through binding and removal of  $Ca^{2+}$  in the surrounding milieu.  $Ca^{2+}$  may regulate several genes responsible for growth and survival of microbes.

## Warnings:

- Do not use the vial if the seal is broken!
- To be used by skilled and trained personnel!
- Dispose unused portions!
- Citra-Lock may not be applied by either direct intravenous injection or added to an infusion!
- Citra-Lock is intended only for use as Catheter Lock!
- Only inject Citra-Lock if the exact catheter volume is known (see instructions catheter manufacturer)!

## Side-Effects:

- No side-effects known, when the solution is used as catheter-lock.
- In the event that Citra-Lock unintentionally passes into the vein, paraesthesia and/or dysgeusia may occur.

## Operating Instructions:

- Aseptic technique must be maintained at all times.
- Prior to initiation of the dialysis session, the Citra-Lock that was instilled in the catheters during the previous session, should be aspirated with the help of a syringe and discarded according to the institution's biohazards waste policy.
- Flush each catheter lumen with 5 ml of sterile 0,9% sodium chloride solution to remove any blood remaining from the previous operation.
- Extract the exact priming volume (prescribed by the catheter manufacturer) of Citra-Lock from the vial, using a 5 ml or smaller syringe (to ensure accurate volume).
- Inject the priming volume of Citra-Lock slowly into the catheter (taking 5-10 sec.).

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Citra-Lock™

Catheter Lock solution

Anti-microbial Anti-clotting Biofilm reduction

Dirinco® B.V. Saffierborch 14 5241 LN Rosmalen The Netherlands [www.citra-lock.com](http://www.citra-lock.com) [info@citra-lock.com](mailto:info@citra-lock.com)

6234/12610052/0211